

Voice Selector Study and wired earphones

Safety Instructions

How to use the device:

Voice Selector Study including a wired earphones is a unique assistive listening device, reducing ambient noise and focusing audio in the selected direction.

Storage and operation temperatures should be between 10° and 45°C, and humidity 10% to 95%.

How to clean the device:

The device should be cleaned with a damp cloth.

Instructions for disposing of electronic equipment:

Electronic waste must be treated according to local regulations. Nuance Hearing Ltd. recommends recycling the device at the regional disposal of electronic equipment. Local laws must be followed.

According to the regulations, we are required to explain the following warnings:

- Warning: Prolonged listening at high volume is not recommended and may cause irreversible damage to your hearing
- Warning: Use of accessories, components or spare parts not supplied by Nuance Hearing may impact use and void the warranty
- The battery in Voice Selector Study is not replaceable and should not be replaced

Standards and regulations

Voice Selector Study meets the strictest international standards for electrical safety and electromagnetic compatibility.

